
U52001 Hebel

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9/04 ALF

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
Instruction sheet

9/04 ALF

- ① Stand base
- ② Stand rod
- ③ Scale
- ④ Pointer
- ⑤ Weights
- ⑥ Holder with ball bearing and knurled nut
- ⑦ Lever arm
- ⑧ Thumb screw

Equipment set to demonstrate the law of the lever, used in experiments on equilibrium and beam-type model scales.

1. Description, technical data

The equipment kit consists of a stand upon which a lever arm is mounted on a ball bearing. There are drilled holes in the middle of the arm running along the arm's axis. These are used to suspend weights from in order to shift the fulcrum's position. A red and white block scale permits a rapid reading of the lever arm's length. A scale with pointer attached to the stand permits precise display of the state of equilibrium.

Lever length:	1 m
Lever weight:	0.458 kg
Number of drill holes:	21
Gap between holes:	50 mm

1.1 Scope of supply

- 1 Lever arm (7)
- 1 Stand rod (2)
- 1 Stand base (1)
- 1 Scale (3)
- 1 Pointer (4)
- 1 Holder with ball bearing and knurled nut (M8) (6)
- 2 Thumb screws (M6)
- 1 Screw
- 10 x 50 g weights (5)

2. Operating principle

The lever belongs to the class of simple machines. The machine we are normally dealing with here is a straight or angled bar. The lever can pivot around an axis (fulcrum) and utilize forces operating at a distance from the fulcrum. Each force results in a torque, which is all the greater in magnitude, the greater the force and the longer the force's line of action is to the rotation axis (length of the arm). A lever with its fulcrum at position O, upon which forces F and G are acting, is in equilibrium, if the left hand torque is equal to the right hand torque (law of levers: force multiplied by the force arm equals load multiplied by the load arm). The following holds true for the forces F , G and the lever arm lengths L_1 , L_2 :

$$G \cdot L_1 = F \cdot L_2$$


Fig. 1: Two arm lever

If the fulcrum of the lever is exactly between the application points of the two forces, this is called a two-arm lever (Fig. 1), if the forces are located on the same side this is termed a single-arm lever (Fig. 2).


Fig. 2: Single-arm lever

3. Operation

3.1 Assembly of the apparatus

- Slide the scale ③ down to the middle of the stand rod ② and fasten it there using the knurled screw.
- Place the stand rod ② onto the stand base ① and use the screw to fasten it.
- Insert the holder with ball bearing ⑥ into the opening of the stand rod and clamp it down tight using the thumb screw ⑧.
- Fasten the lever arm ⑦ to the ball bearing using the knurled nut ⑥.
- Screw the pointer ④ into the knurled nut ⑥.
- Hang the weights at the desired locations.

4. Sample experiments

4.1 Verify the law of levers using the two-arm lever

- Fasten the lever arm in the middle of the ball bearing.

- Hang 5 50 g weights (each approx. 2.5 N) on the left side of the lever arm in the third bore hole from the center.
- Hang 3 50 g weights (each approx. 1.5 N) on the right side of the lever arm in the fifth bore hole from the center.
- The lever is in a state of equilibrium.
- Repeat the experiment using other weight combinations.

4.1 Verify the law of levers using a single-arm lever

Additionally required equipment:


One dynamometer 2 N, 5 N and 10 N each (e.g. U20033/U20034/U20035)

- Fasten the lever arm to the ball bearing in the last bore hole (see Fig. 2).
- Hang 5 50 g weights (each approx. 2.5 N) on the right side of the lever arm at the location of force G .
- Instead of force F deploy a dynamometer.
- Bring the lever into a state of equilibrium and then note down the force registered at the dynamometer.
- In this experiment the weight of the lever arm must be taken into consideration (0.458 kg corresponds to 4.49 N).

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